

THE NON-CONTACT OPERATION PREVENTS CONTAMINATION



CAT No. I-WMSO-2

INSTRUMED's **Wall Mounted Scrub Station**

**Automatic Aseptic Stainless Steel Scrub-up
For Operation Theatres, Cathlabs.
& Allied Sterile Areas**



INSTRUMED (India) International

**WALL
MOUNTED
SAFE
HYGIENIC
CONVENIENT
ELEGANT**



**NO LEAKAGE
& DRIPPING
OF TAPS
DIFFERENT
MODELS &
OPTIONS TO
CHOOSE**

INSTRUMED's STAINLESS STEEL AUTOMATIC ASEPTIC SCRUB STATION

CONSTRUCTION: The entire unit is fabricated of high quality stainless steel sheets & pipes. The design of trough ensures no stagnation of water. Maintenance door/s in front. A manual mixer is provided for mixing Hot & Cold water. The unit is compatible with all standard & regular plumbing or can be fabricated as per site conditions.

OPERATION: Totally hand free operation. The flow of water is controlled by Sensors.

EMERGENCY BACKUP FACILITY: Elbow operated tap is provided which can be used in case of Power/ Sensor failure.

OPTIONAL ADDITIONAL FACILITY: (In Sensor Operated Models) An additional high speed "JET FLOW" can be provided for washing of rigid and flexible endoscopes, instruments, tubes etc.

SITE CONDITION REQUIREMENTS: Clean site with water inlet & outlet (Drain) & electrical points.

OPTIONAL ADDITIONAL FACILITY : Another facility or feature required can be incorporated - if possible.

MODELS:

(A) SENSOR OPERATED:

CAT NO. I-WMSO-1 for Single Person
I-WMSO-2 for Two Persons To Scrub Simultaneously
I-WMSO-3 for Three Persons To Scrub Simultaneously

(B) MANUALLY OPERATED: (The flow of water is regulated by Elbow Operated Tap/s.).

CAT NO. I-WMM-1 for Single Person
I-WMM-2 for Two Persons To Scrub Simultaneously
I-WMM-3 for Three Persons To Scrub Simultaneously

WARRANTY: All units carry a warranty of 12 months against any manufacturing defect.

NOTE: We reserve full & absolute rights to change/alter the Models/Specification for improvement of the system without notice.



INSTRUMED (India) International

An ISO 9001-2000 Certified Co.

4852/24, Ansari Road, Darya Ganj, New Delhi- 110 002 (India).

Ph : (91)-011-23277175, 23256681 Telefax : (91)-011-23256681

E-mail : info@instrumedindia.com Website : www.instrumedindia.com